

LSINC1114 – Analysis of biological data

Session 8 – Medical imaging and DICOM


Sébastien Jodogne

EPL / ICTEAM
UCLouvain, Belgium

Academic year 2022-2023

1. Brief overview of medical imaging
2. DICOM
3. HTML5 canvas

Brief overview of medical imaging



- Belgium (2013): 33 millions of imaging studies for 11 millions of persons.
- Reason: **Multimodal and longitudinal** imaging (oncology, cardiovascular diseases, surgery, neurology...)

- Many lesions/disorders are **internal** → How can we detect them and select the best treatment?
- How can we see **inside the body in a non-invasive way**, i.e. without “opening” it with a scalpel?



Invasive



Non-invasive

Credits: B. Macq, J. Lee, F. Peeters, G. Kerckhofs (LGBIO2050 @ UCLouvain)


It all began in **1895** with Wilhelm Conrad **Röntgen**'s picture of his wife Bertha's hand:



Credits: B. Macq, J. Lee, F. Peeters, G. Kerckhofs (LGBIO2050 @ UCLouvain)

How does this work?

- An homogeneous beam of X-rays is produced by an **X-ray generator** and is directed towards an object.
- According to the density and composition of the different components of the object, **X-rays are more or less absorbed**.
- X-rays are then captured behind the object by a **detector**, which gives a 2D image as if it was transparent somehow.



Credits: B. Macq, J. Lee, F. Peeters, G. Kerckhofs (LGBIO2050 © UCLouvain)

- Conventional X-ray photographic plate, to be chemically developed, then scanned (**RadioGraphic imaging = RG**).
- Fluorescent screen, recorded by a visible-light camera (**Radio Fluoroscopy = RF**).
- Better: Photo-stimulable phosphor plate, where excited electrons are read using a laser, which is reusable (**Computed Radiography = CR**).
- Semiconductor detectors can directly create an image, since 2000's (**Digital Radiography = DR**).



RF





CR



DR




Radiography device rotating around the patient
→ **Tomographic reconstruction** (inverse Radon transform, 1917) → 3D image



1971: EMI scanner
(Sir Godfrey Hounsfield)


<https://thelongandshort.org/enterprise/how-the-beatles-revolutionised-medical-imaging>




Ring artifact




Ring artifact



Streak artifact (metal)




10 mA/10 ms





80 mA/10 ms

Noise / partial volume effect (low-dose)




Motion artifact

- Other waves and/or techniques are used in medical imaging:
 - DR - Digital radiography → X-ray transmission
 - CT - Computed tomography → X-ray transmission
 - CBCT - Cone beam computed tomography → X-ray transmission
 - MRI - Magnetic resonance imaging → magnetic fields
 - US - Echography → ultrasound waves
 - SPECT - Single photon emission computed tomography → gamma ray emission
 - PET - Positron emission tomography → gamma ray emission
- Why multiple modalities?
 - Different purposes: anatomical vs. functional.
 - Different properties: resolution, contrast, etc.



Credits: B. Macq, J. Lee, F. Peeters, G. Kerckhofs (LGBIO2050 @ UCLouvain)


Electromagnetic spectrum



- Wavelength (in vacuum): $\lambda = c/f$, where:
 - $c = 299,792,458[m/s]$ is the speed of light in vacuum.
 - f is the frequency of the wave [Hz] = [$1/s$].
 - Unit: meter.
- Photon energy: $E = h \cdot f$, where:
 - $h = 6.626 \cdot 10^{-34}[J \cdot s]$ is Planck's constant.
 - Unit: Joule [J].
 - Electronvolt [eV] = $1,6 \cdot 10^{-19}[J]$: Amount of kinetic energy gained by a single electron accelerating from rest through an electric potential difference of one volt in vacuum.
 - Photon energy expressed in electronvolts: $E \approx 4,1 \cdot 10^{-15} \cdot f$.
- *Don't learn those figures ;)*


https://en.wikipedia.org/wiki/Electromagnetic_spectrum

DICOM




“One **OPEN standard** to rule them all...”

“Big imaging” in hospitals



Everything is driven by software
using the DICOM open standard!

Also “small imaging” and outside of hospitals!



DICOM is everywhere!



File format

```
<?xml version="1.0"?>
<person id="001">
  <name>Kris</name>
  <address>
    <street>...</street>
    <city>...</city>
  </address>
</person>
```

Acquisition data
(kind of XML)





One slice of image
(kind of PNG)

Network protocol



Early example of Web service:
« store / query / retrieve »



Large variations in the volume of one medical image: radiography – 10MB, mammography – 100MB, CT-scan – 500MB, histology – 10GB


- “Digital Imaging and Communications in Medicine”.
- De-facto **open** standard for digital medical imaging in all hospitals around the world.
- First normalized in **1985** by the ACR (*American College of Radiology*) and NEMA (*National Electrical Manufacturers Association*).
- Goal: Maximal technical **interoperability** between hardware and software, independently of their vendors.
- This session is about the DICOM file format.

Libraries



dcm4che.org

Viewers



Starviewer
Medical Imaging Software

Servers


dcm4chee



conQEST


ORT[H]ANC

digital
wallonia
.be



Ancillary PACS for research or automation (yet main PACS)

Inter-site exchanges (including to the cloud)



Teleradiology portals

- A DICOM file is made of three parts:
 - The “**meta-header**” that identifies the file,
 - The “**dataset**” that contains the acquisition data,
 - The “**pixel data**” that contains the actual 2D image.
- The pixel data can be **compressed** according to a “*transfer syntax*”:
 - In the clinical workflow, **no compression (raw data)** for maximum compatibility.
 - **Lossless** (non-destructive) compression:
 - RLE, JPEG2k, JPEG-LS.
 - Default choice for medical traceability.
 - **Destructive** compression:
 - JPEG.
 - Used for teleradiology (low-resolution image to reduce network bandwidth), and possibly digital pathology or visible-light imaging.
- A DICOM file can be **multi-frame**:
 - Several 2D images are concatenated inside the pixel data.
 - Commonly used by US devices (cine), PET-scanners and radiotherapy doses.

<https://book.orthanc-server.com/dicom-guide.html#dicom-file-format>

- The meta-header and the dataset are lists of **DICOM tags** associated with values:
 - A DICOM tag is a **pair of hexadecimal numbers** (up to 65535): The **group** and the **element**.
 - The value associated with a DICOM tag has a data type: the **value representation**.
 - The DICOM standard lists mandatory (“types 1 and 2”) or optional (“type 3”) DICOM tags for each imaging modality.
 - DICOM tags can be nested (recursive structure like HTML or JSON) → “**sequences**”.
- The DICOM standard defines a **nomenclature of DICOM tags**, for instance:
 - (0x0010,0x0010) → Patient's name
 - (0x0010,0x0030) → Patient's birth date
 - (0x0008,0x1030) → Study description
 - (0x0008,0x1010) → Station name
 - (0x0008,0x0060) → Modality (MR, CT, RX, US...)
 - (0x0028,0x0010) and (0x0028,0x0011) → Rows and columns (height/width)
- The **pixel data** is actually part of the dataset, as the (0x7fe0,0x0010) tag.

<https://book.orthanc-server.com/dicom-guide.html#dicom-file-format>

NB: The command-line tools “dcm2xml” (from DCMTK) and “gdcmxml” (from GDCM) can also be used.

Extremely simple: Concatenate row by row, then frame by frame!

0,0	0,1	0,2
1,0	1,1	1,2
2,0	2,1	2,2



0,0	0,1	0,2	1,0	1,1	1,2	2,0	2,1	2,2
-----	-----	-----	-----	-----	-----	-----	-----	-----

Image of size “3x3” (single-frame)






Image of size “3x3x3” (multi-frame)





Large variations in the volume of one imaging study:

- radiography ~ 10MB (2 files of 5MB),
- mammography ~ 100MB (4 files of 30MB),
- CT-scan ~ 500MB (1000 files of 500KB),
- MRI ~ 1GB (large variations).

- 3D volumes are typically encoded as multiple DICOM instances (one per 2D slice).









- “Russian dolls” model.
- “Studies” vs. “series”:
 - Nuclear medicine → two series (CT + PET) that are part of the same study.
 - CT-scanner → multiple series (scout, sagittal, coronal, axial, preview, dose reports...).
 - MRI → multiple series depending on the acquisition sequence (T1, T2, FLAIR, BOLD...).
- DICOM specifies one **module** (set of DICOM tags) for each of the “Patient”, “Study”, “Series” and “Instance” levels.

DICOM hierarchy, as displayed by Orthanc





- Each DICOM instance contains 4 **identifier tags**:
 - (0x0010,0x0020) → Patient ID (in the “Patient” module)
 - (0x0020,0x000d) → Study Instance UID (in the “Study” module)
 - (0x0020,0x000e) → Series Instance UID (in the “Series” module)
 - (0x0008,0x1018) → SOP Instance UID (in the “Instance” module)
- Those can be used as **primary keys** for the database of a PACS.
- Uniqueness:
 - Study/Series/SOP instance UID are guaranteed to be **globally unique** (i.e. all around the world).
 - Patient ID is only unique **locally unique** (inside the hospital)
→ need for **reconciliation** if medical images are transferred from one hospital to another.



Table A.3-1. CT Image IOD Modules

IE	Module	Reference	Usage
Patient	Patient	C.7.1.1	M
	Clinical Trial Subject	C.7.1.3	U
Study	General Study	C.7.2.1	M
	Patient Study	C.7.2.2	U
	Clinical Trial Study	C.7.2.3	U
Series	General Series	C.7.3.1	M
	Clinical Trial Series	C.7.3.2	U
Frame of Reference	Frame of Reference	C.7.4.1	M
Equipment	General Equipment	C.7.5.1	M
Image	General Image	C.7.6.1	M
	Image Plane	C.7.6.2	M
	Image Pixel	C.7.6.3	M
	Contrast/Bolus	C.7.6.4	C - Required if contrast media was used in this image
	Device	C.7.6.12	U
	Specimen	C.7.6.22	U
	CT Image	C.8.2.1	M
	Overlay Plane	C.9.2	U
	VOI LUT	C.11.2	U
	SOP Common	C.12.1	M
	Common Instance Reference	C.12.2	U

Table C.7-3. General Study Module Attributes

Attribute Name	Tag	Type	Attribute Description
Study Instance UID	(0020,000D)	1	Unique identifier for the Study.
Study Date	(0008,0020)	2	Date the Study started.
Study Time	(0008,0030)	2	Time the Study started.
Referring Physician's Name	(0008,0090)	2	Name of the patient's referring physician
Study ID	(0020,0010)	2	User or equipment generated Study identifier.
Accession Number	(0008,0050)	2	A RIS generated number that identifies the order for the Study.

- DICOM standard groups DICOM tags by **modules**.
 - Each type of imaging modality is associated with a **IOD** that lists the mandatory and optional modules.
 - One module can be shared by multiple IOD.
-
- A module can be mandatory (M), optional (U) or conditional (C) in an IOD.



- Very large 2D images → $100,000 \times 100,000$ pixels → up to 20GB.
- Different zoom levels are precomputed (**multi-scale pyramid**).
- Stored as many multi-frame files, each tile being one frame.
 - A single tile will typically be 256×256 or 512×512 pixels.

DICOM encoding
(cf. Google Maps)

<https://wsi.orthanc-server.com/demo/>

```
import java.io.FileInputStream;
import java.io.IOException;
import org.dcm4che3.data.Attributes;
import org.dcm4che3.data.Tag;
import org.dcm4che3.io.DicomInputStream;

public class Test {
    public static void main(String[] args) throws IOException {
        FileInputStream stream = new FileInputStream("hand.dcm");

        Attributes dataset;

        try (DicomInputStream din = new DicomInputStream(stream)) {
            dataset = din.readDataset();
        }

        System.out.println("Patient name: " + dataset.getString(Tag.PatientName));
        System.out.println("Study date: " + dataset.getString(Tag.StudyDate));
    }
}
```

- Java applications will typically use the **dcm4che** library.
- Class `be.uclouvain.DicomImage` provides an easy-to-use wrapper.

HTML5 canvas

From HTML:

```
<canvas id="my-canvas" width="640" height="480"></canvas>
```



Our BEST Rendering Library has to wrap `<canvas>` within a `<div>` to control its size using CSS:

```
<div id="my-drawing" class="container">
  <canvas></canvas>
</div>
```

From JavaScript:

```
var canvas = document.createElement('canvas');
canvas.width = 640;
canvas.height = 480;
```

- This creates an invisible “offscreen” canvas
- `canvas` is typically a global variable


Drawing on a canvas is
done by getting its
rendering context:

```
// HTML canvas
var ctx = document.getElementById('my-canvas').getContext('2d');

// Offscreen canvas
var canvas = document.createElement('canvas');
canvas.width = 640;
canvas.height = 480;
var ctx = canvas.getContext('2d');
```

<https://developer.mozilla.org/en-US/docs/Web/API/CanvasRenderingContext2D>

```
// Set line width  
ctx.lineWidth = 10;  
  
// Wall  
ctx.strokeRect(75, 140, 150, 110);  
  
// Door  
ctx.fillRect(130, 190, 40, 60);  
  
// Roof  
ctx.beginPath();  
ctx.moveTo(50, 140);  
ctx.lineTo(150, 60);  
ctx.lineTo(250, 140);  
ctx.closePath();  
ctx.stroke();
```



<https://developer.mozilla.org/en-US/docs/Web/API/CanvasRenderingContext2D>

	0	1	2	3
0	23	216	120	55
1	4	89	158	130
2	65	76	189	34
3	19	234	7	45

x
y

	0	1	2	3
0	23	216	120	55
1	4	89	158	130
2	65	76	189	34
3	19	234	7	45

c
r

	1	2	3	4
1	23	216	120	55
2	4	89	158	130
3	65	76	189	34
4	19	234	7	45

c
r

	y	3	2	1	0
0		23	216	120	55
1		4	89	158	130
2		65	76	189	34
3		19	234	7	45

x

DICOM / HTML5 /
BufferedImage in Java

numpy in Python /
RealMatrix in Java

Matlab / Octave

Traditional
mathematical
coordinates

<https://doi.org/10.1007/978-3-030-39364-9> (Section 2.4.1)

```
// Change color of stroke (lines)
ctx.strokeStyle = 'red';

// Change color of filling
ctx.fillStyle = 'gray';

// Drawing a circle
ctx.beginPath();
ctx.arc(centerX, centerY, 100, 0, 2.0 * Math.PI);
ctx.stroke();

// Filling a circle
ctx.beginPath();
ctx.arc(centerX, centerY, 100, 0, 2.0 * Math.PI);
ctx.fill();
```

<https://developer.mozilla.org/en-US/docs/Web/API/CanvasRenderingContext2D>

From an offscreen canvas:

```
<canvas id="my-canvas" width="640" height="480"></canvas>
```



```
var canvas = document.createElement('canvas');
// ...Draw onto the offscreen canvas...

var ctx = document.getElementById('my-canvas').getContext('2d');
ctx.drawImage(canvas, x, y);
```

From a PNG/JPEG “sprite” loaded by HTML:



```
<canvas id="my-canvas" width="640" height="480"></canvas>

```




```
var sprite = document.getElementById('sprite');
var ctx = document.getElementById('my-canvas').getContext('2d');
ctx.drawImage(sprite, x, y);
```

From image to canvas coordinates




t_x = shift in x
 t_y = shift in y



$$\begin{cases} x' = x + t_x \\ y' = y + t_y \end{cases}$$

s_x = scaling factor in x

s_y = scaling factor in y




- NB: Before applying zoom, the center of the drawing should generally be translated at coordinate $(0, 0)$.

$$\begin{cases} x' = a_{11} \cdot x + a_{12} \cdot y + a_{13} \\ y' = a_{21} \cdot x + a_{22} \cdot y + a_{23} \end{cases} \Leftrightarrow \begin{pmatrix} x' \\ y' \\ 1 \end{pmatrix} = \underbrace{\begin{bmatrix} a_{11} & a_{12} & a_{13} \\ a_{21} & a_{22} & a_{23} \\ 0 & 0 & 1 \end{bmatrix}}_{A \in \mathbb{R}^{3 \times 3}} \begin{pmatrix} x \\ y \\ 1 \end{pmatrix}$$

↑ "2D homogeneous coordinates"


+



rotation shearing

zoom translation


Chaining transformations



Matrix multiplication


$$\begin{pmatrix} x' \\ y' \\ 1 \end{pmatrix} = A \begin{pmatrix} x \\ y \\ 1 \end{pmatrix}; \begin{pmatrix} x'' \\ y'' \\ 1 \end{pmatrix} = A' \begin{pmatrix} x' \\ y' \\ 1 \end{pmatrix}; \begin{pmatrix} x''' \\ y''' \\ 1 \end{pmatrix} = A'' \begin{pmatrix} x'' \\ y'' \\ 1 \end{pmatrix}$$
$$\Rightarrow \begin{pmatrix} x''' \\ y''' \\ 1 \end{pmatrix} = \boxed{A'' A' A} \begin{pmatrix} x \\ y \\ 1 \end{pmatrix}$$

Reversing a transformation



Matrix inversion

$$\begin{pmatrix} x' \\ y' \\ 1 \end{pmatrix} = A \begin{pmatrix} x \\ y \\ 1 \end{pmatrix} \Rightarrow \begin{pmatrix} x \\ y \\ 1 \end{pmatrix} = \boxed{A^{-1}} \begin{pmatrix} x' \\ y' \\ 1 \end{pmatrix}$$



Chapter 10

Auxiliary function:

```
var CreateTransform = function(a11, a12, a13,
                               a21, a22, a23) {
    // "DOMMatrix" uses a transposed matrix
    return new DOMMatrix([ a11, a21,
                          a12, a22,
                          a13, a23 ]);
};
```

```
var ctx = document.getElementById('my-canvas').getContext('2d');
var backup = ctx.getTransform();

var zoom = CreateTransform(2, 0, 0,
                           0, 2, 0);
ctx.setTransform(zoom);

var pan = CreateTransform(1, 0, 100,
                        0, 1, 200);
ctx.setTransform(pan);





// This is the matrix product "zoom * pan"
var firstPanThenZoom = zoom.multiply(pan);
ctx.setTransform(firstPanThenZoom);


ctx.setTransform(backup); // Back to the original transformation
```



JavaScript library developed as a support to the BEST summer course
“Hacking medical imaging” in Louvain-la-Neuve (July 2022).

Mouse conventions in radiology







Scene coordinates

The designer of a Web application for medical imaging wants to focus on drawing the **content of scene**, not on the coordinates of the HTML page!


First wish: Fit scenes to browser window






Scene coordinates

Viewport coordinates (zoom/pan)



Browser coordinates
(size of the window)

Chain of 3 affine transformations




The **BEST** rendering library
manages this for you!



NB: Don't hesitate to look at the source code "`BestRendering.js`"
to know how transformations and mouse events are handled!

Example (1/4): The “magic” rectangle

Scene coordinates



Example (2/4): Drawing a simple scene

index.html

```
<!DOCTYPE html>
<html>
  <head>
    <meta charset="utf-8">
    <style>
      html, body {
        width: 100%;
        height: 100%;
        margin: 0px;
        border: 0;
        overflow: hidden; /* Disable scrollbars */
        display: block; /* No floating content on sides */
      }
      #my-drawing {
        position: absolute;
        left: 0;
        top: 0;
        width: 100%;
        height: 100%;
      }
    </style>
  </head>
  <body>
    <div id="my-drawing" class="container">
      <canvas></canvas>
    </div>
    <script src="BestRendering.js"></script>
    <script src="app.js"></script>
  </body>
</html>
```

app.js

```
var color = 'red';

function GetExtent() {
  // Define the extent of the scene, i.e. (x1,y1) and (x2,y2)
  return BestRendering.CreateExtent(-1, -1, 2, 1);
}

function Draw(ctx) {
  // "ctx" is now working in scene coordinates!
  ctx.fillStyle = color;
  ctx.fillRect(-1, -1, 3 /* width */, 2 /* height */);


  ctx.strokeStyle = 'black';
  ctx.moveTo(0.5, -1); ctx.lineTo(0.5, 1); ctx.stroke();
  ctx.moveTo(-1, 0); ctx.lineTo(2, 0); ctx.stroke();
}

BestRendering.InitializeContainer('my-drawing', Draw, GetExtent);
```




Example (3/4): Re-fit the canvas after zoom/pan

```
document.addEventListener('keydown', function(event) {  
    if (event.key == 'f') {  
        BestRendering.FitContainer('my-drawing');  
    }  
});
```



Example (4/4): Redraw after a change in the scene

```
document.addEventListener('keydown', function(event) {
    if (event.key == 'g') {
        color = 'green';
        BestRendering.DrawContainer('my-drawing');
    } else if (event.key == 'r') {
        color = 'red';
        BestRendering.DrawContainer('my-drawing');
    }
});
```



“g”



“r”



